



## RAPID TOX CUP II™

- **Easy to use.** No timing required: control line indicates test is ready to be interpreted
- **Fast:** results in 5 minutes
- **Accurate:** >99% correlation to GC/MS at 95% confidence level, means: Results you can trust.
- **Reliable:** Tested against over 600 compounds for cross reactivity
- Can detect up to 16 drugs
- Customization of test is available to meet your specific drug testing needs. The following strip assays (in combination or alone) are available to provide standard configurations of 4, 8 or 12 panels that detect:
  - Strip 1: OPI 300, METH 500, COC 150, THC
  - Strip 2: BAR, BZO, PCP, AMP 500
  - Strip 3: MDMA 500, MTD, OXY, BUP
  - Can add: SVT, TRAM, FENT, ETG, K2 and more
- **Call or visit our website today for more information!**

*For results you can count on, count on us.*

ABMC is an ISO certified manufacturing company that manufactures accurate, cost effective, easy to use, onsite drug of abuse tests. Our clients use our tests with confidence because of our proven accuracy and longevity in the industry. With an extensive portfolio of urine and oral fluid drug screens providing results in minutes, one can test for a combination of 1-14 different drugs with custom configurations available. Web based or onsite training, with certification, is available upon request. Technical Support is available 24/7/365.



ONLINE & IN PERSON TRAINING



MADE IN USA



24/7/365 TECHNICAL SERVICE AVAILABLE



DOCUMENTED ACCURACY

122 SMITH ROAD

KINDERHOOK, NEW YORK 12106

WWW.ABMC.COM

800.227.1243

518.758.8158