



## CONVENIENCE

- Easy to use: integrated collection and testing
- Temperature strip on cup
- No timing required: control line indicates read time
- Results ready in 3-5 minutes
- Collection cup allows for visual integrity check
- Closed specimen cup: limits exposure to specimen
- Ship in same container for confirmation

## PERFORMANCE

- Built-in procedural control
- Results stable for up to one hour
- >99% correlation to GC/MS at 95% confidence level
- Colloidal gold technology
  - = superior sensitivity, accuracy & precision
- Tested for cross-reactivity against over 600 compounds

SAMHSA cut-off levels

Integrated test strips are US FDA 510(k) cleared

## FLEXIBILITY

- 2- 12 drugs of abuse per cup
- 15 drug assays to choose from
- Specimen validity test available
- Customized configurations available

Multiple cut-off options for certain drugs

## TRAINING

- Web based training with certification
- Onsite training

## TECHNICAL SUPPORT

Available 24/7/365



**AMP**  
amphetamines  
1000ng/mL

**BAR**  
barbiturates  
300ng/mL

**BZO**  
benzodiazepines  
300ng/mL

**BUP**  
buprenorphine  
12.5ng/mL

**COC**  
cocaine  
150 & 300ng/mL

**MDMA**  
ecstasy  
1000ng/mL

**MTD**  
methadone  
300ng/mL

**METH**  
methamphetamines  
1000ng/mL

**NIC\***  
nicotine  
200ng/mL

**OPI**  
opiates  
300 & 2000ng/mL

**OXY**  
oxycodone  
100ng/mL

**PCP**  
phencyclidine  
25ng/mL

**PPX**  
propoxyphene  
300ng/mL

**THC**  
marijuana  
50ng/mL

**TCA**  
tricyclic  
antidepressants  
1000ng/mL

