

RAPID STAT™



CONVENIENCE

Easy to use

All collections can be observed

No restroom needed

No shy bladder issues

Specimen collected in just minutes

No specimen dilution or adulteration

Volume indicator: ensures correct specimen volume

Rapid results=quick decisions

Designed to make confirmation testing simple

PERFORMANCE

Built-in procedural control

Colloidal gold technology=
superior sensitivity accuracy and precision

Results in about 10 minutes

Tested for cross-reactivity against most common
drugs and other substances

Detects parent THC

SAMHSA states: "The only way to detect marijuana use is through the presence of the parent drug (THC) in the oral fluid....."

FLEXIBILITY

Eight drugs available

Multiple drugs tested simultaneously

Multiple drug combinations available

TRAINING

Web based training with certification

Onsite training

TECHNICAL SUPPORT

Available 24/7/365

AMP
amphetamines
25ng/mL

BZO
benzodiazepines
10ng/mL

COC
cocaine
12ng/mL

MTD
methadone
15ng/mL

METH/MDMA
methamphetamine
25ng/mL

OPI
opiates
10ng/mL

PCP
phencyclidine
2.5ng/mL

THC
marijuana
25ng/mL