



## CONVENIENCE

Easy to use: integrated collection and testing

Temperature strip on cup

No timing required: control line indicates read time

Results ready in 3-5 minutes

Collection cup allows for visual integrity check

Closed specimen cup: limits exposure to specimen

Ship in same container for confirmation

## PERFORMANCE

Built-in procedural control

Results stable for up to 6 hours

>99% correlation to GC/MS at 95% confidence level

Colloidal gold technology

= superior sensitivity, accuracy & precision

Tested for cross-reactivity against over 600 compounds

SAMHSA cut-off levels

US FDA 510(k) Cleared

## FLEXIBILITY

Up to 15 drugs of abuse per cup

15 drug assays to choose from

Specimen validity test available

Customized configurations available

Multiple cut-off options for certain drugs


## TRAINING

Web based training with certification

Onsite training

## TECHNICAL SUPPORT

Available 24/7/365



**AMP**  
amphetamines  
500 & 1000ng/mL

**BAR**  
barbiturates  
300ng/mL

**BZO**  
benzodiazepines  
300ng/mL

**BUP**  
buprenorphine  
12.5ng/mL

**COC**  
cocaine  
150 & 300ng/mL

**MDMA**  
ecstasy  
500 & 1000ng/mL

**MTD**  
methadone  
300ng/mL

**METH**  
methamphetamines  
500 & 1000ng/mL

**NIC\***  
nicotine  
200ng/mL

**OPI**  
opiates  
300 & 2000ng/mL

**OXY**  
oxycodone  
100ng/mL

**PCP**  
phencyclidine  
25ng/mL

**PPX**  
propoxyphene  
300ng/mL

**THC**  
marijuana  
50ng/mL

**TCA**  
tricyclic  
antidepressants  
1000ng/mL

\*Nicotine is not  
USA 510(k) cleared

1 2 2 S M I T H R O A D , K I N D E R H O O K , N E W Y O R K 1 2 1 0 6  
W W W . A B M C . C O M ■ 8 0 0 . 2 2 7 . 1 2 4 3 ■ 5 1 8 . 7 5 8 . 8 1 5 8