

The On-The-Spot Leader in Saliva Alcohol Testing



Q.E.D.[®]
SALIVA ALCOHOL TEST

The Q.E.D. Saliva Alcohol Test is the first FDA-cleared, CLIA-waived and DOT-approved point-of-care test for saliva alcohol.

Simple-to-Use

- 3 easy-step procedure
- Results in 2 minutes
- Built-in quality control

Proven Accuracy

- >98% correlation to blood gas chromatography
- Precise quantitative results (0-145 mg/dL)

Trusted and Reliable

- FDA-cleared; CLIA-waived
- DOT-approved
- Accepted by major insurance carriers (CPT 80305)

The Value of a Reliable and Trusted Brand



OraSure Technologies



Simple-to-use


3 EASY STEPS



1. Swab mouth.




2. Insert collector.



3. Read in 2 minutes.

Built-in quality control for correct procedural assurance

Proven clinical lab accuracy at point-of-care



The Q.E.D.® test can be used as an accurate screen in place of blood tests. In clinical trials, saliva alcohol levels measured with the Q.E.D. tests demonstrated high correlation to blood analyzed by gas chromatography ($r=0.98$)¹ and produced no false positive or false negative results².

Data on file. ¹Study conducted by Thomas Jefferson University Hospital, Department of Emergency Medicine. ²Study conducted by Rutgers University, Center of Alcohol Studies.

The Q.E.D. test is used for . . .

DOT and non-DOT worksite & occupational health testing

- Random testing
- Post-accident
- Reasonable cause
- Follow-up testing

Physician offices/emergency & trauma depts.

- Risk-reduction testing
- Pre-discharge
- Admission screening
- Seizure or unconscious patients
- Compliance testing

Substance abuse and criminal justice

- Testing upon admissions
- Interventions
- Counseling sessions
- Random testing

Reimbursed by major insurance carriers

2017 Medicare National Limit Code 80305 QW*	Estimated Level of Reimbursement	
	30 Count Kit	100 Count Kit
\$14.96	\$448.80	\$1,496.00

* The QW modifier is required for Medicare carriers effective 4/1/17 and is retroactive for dates of service beginning 1/1/17. Commercial payors may not require QW. It is the provider's responsibility to review each payor policy for appropriate billing guidelines.



Bethlehem, PA 18015 USA
Made in the USA

Order Information
31150B
31150C
31150S
31050S

Description
Box of 30 tests
Box of 100 tests
DOT Controls
A150 Controls



American Bio Medica Corporation
122 Smith Road
Kinderhook, NY 12106
800-227-1243 / 581-758-8158
www.abmc.com